

# CA-800 Corneal Analyser





### Ease of use

The CA-800 is extremely easy to handle and use. From image acquisition to analysis, the on-board software is intuitive and user-friendly and the 10.1-inch capacitive touch screen provides quick navigation. Visual guidance supports fast and easy alignment and focusing on the eye; the "best image" selection mode automatically acquires the best-focused image. The CA-800 is a placido-based topography system that delivers accurate, high resolution images of the anterior corneal surface. The keratoscope cone with 24 rings equally spaced on a 43D sphere analyses over 100.000 data points, with axial and instantaneous curvature evaluation.

#### **Integrated PC**

The brand new compact design of the CA-800 includes a fully integrated PC, so that an external PC is not required to manage a patient database for archiving and re-analysis purposes. The patient database is stored on an internal 320Gb SATA hard disk and the CA-800 includes a 32Gb SSD for a quick startup of the instrument and user interface.

# Accurate, full examination of the anterior corneal surface



## CA-800 fully featured

#### » Topography map

- » Map full screen mode
- » Ring editing
- » Keratoconus screening (KPI)
- » Full 3D map of corneal surface
- » Automated best image selection
- » OD/OS results on same screen

#### » Corneal wavefront (Zernike) analysis

- » Corneal surface height map
- » Comparison map
- » Reviewing of previous patient examinations
- » Differential map
- » Post-operative monitoring of corneal healing

#### » Pupillometry

- » Automated pupil recognition
- » Dynamic, Photopic, Mesopic & Scotopic
- » Latency graph
- » Real time fluorescein acquisition and imaging
   » Internal yellow barrier filter
- » White to white measurement
- » Meibomian gland analysis

#### » Contact lens fitting simulation

- » Complete contact lens fitting software
- » Contact lens database on-board


#### » Toric IOL calculation

- » Oculentis
- » 19" Capacitive touch screen
- » Fully integrated patient database







## All features accessible on just one screen












Acquisition The CA-800 is easy to use. Visual signals support fast and easy alignment and focusing on the patients eye. The CA-800 has a Right and Left eye detection and prevents incorrect savings in Right / Left eye measurements. The automated best image selection mode in the software of the CA-800 decides the best focused position and automatically acquires the image. Acquistions can be made for topography, pupillometry and real time fluorescein imaging.

#### **Keratoconus screening**

With the CA-800, signs of asymmetry of the cornea can easily be detected even in an early stage. By analyzing the apical curvature, apical gradient and symmetry of the cornea, a Keratoconus probability index will show in colour code (green, yellow & red) if the topography is compatible with Keratoconus. With the CLMI (Cone Location and Magnitude Index) it is easy to follow-up on keratoconus and keratoconus-like patterns.




The Zernike analysis module consists of 36 polynomials into the 7th order, and provides a clear view on the optical deficiencies which can disturb vision. Based on this information, the CA-800 provides the visual acuity summary. Zernike analysis is the basis for the calculation of the ablation area for laser treatment. The Zernike expansion coefficient is used to determine which component(s) dominate the aberration structure of the cornea and to what degree.

# CA-800 - Corneal Analyser



#### **Corneal Zernike analysis**



#### **Corneal comparison & differential map**

With the CA-800, it is easy to compare topography maps between two examinations of a patient, which can be used for follow up and for pre- and post-operative corneal analysis. With the differential map, progress in recovery of the cornea can be observed after refractive surgery. Parameters such as keratometry, apical curvature and corneal symmetry can be analysed to follow the development of any corneal surface changes. The CA-800 comparison and differential maps help you with the treatment of collagen cross linking to stop the development of corneal keratoconus.

# Pupillometry

The CA-800 is equipped with two white LED's for dynamic and static pupillometry. With the CA-800 on-board, the user can check the pupil position and diameter (from Photopic to Scotopic condition) in relation to the position of the optical zone in Ortho-K, laser treatment or refractive surgery treatments. Dynamic pupillometry provides clear information on the reaction time of the pupil and the contraction of the pupil.



### Meibomian gland analysis

With the Infra-red illumination of the CA-800, the Meibomian Glands of the upper and lower eyelid can be captured and analysed. Posterior blepharitis is the most common form of lid margin disease. MGD (Meibomian Gland Dysfunction) can cause or exacerbate dry eye symptoms and eyelid inflammation. The oil glands become blocked with thickened secretions. Chronically clogged glands eventually become unable to secrete oil which results in permanent changes in the tear film and dry eyes. With the CA-800, MGD can easily be observed and compared with previous Meibomian gland examinations of the patient.



#### Fluorometry

The CA-800 incorporates eight blue LED's for fluorescein images and real time fluorescein videos which are essential for contact lens fitting. During every measurement, the CA-800 automatically registers the pupil diameter, which is critical information during contact lens fitting. Real time fluorescein films allow the eye care practitioner to judge the movement of the contact lens on the cornea, the distribution of the tear film under the contact lens as well as the wetting of the outer contact lens surface. The corneal condition can be observed by recording a real time fluo film without wearing a contact lens. The tear film condition, corneal artifacts and break up tear time (BUT) can be observed.




# CA-800 - Corneal Analyser



#### **Contact lens fitting simulation**

The CA-800 provides the perfect platform for contact lens fitting. Simulation software is provided on-board, which automatically selects the best fitting contact lens based upon an included complete contact lens database for all the main manufacturers (upgradable and customizable by the user). With the option to input refractive powers, the contact lens proposal is accurate and complete. The on-board fluorescein acquisition system allows full control of the contact lens position on the eye. The comparison between different contact lenses is easy in order to ensure the best fit.


# Reports Topography




ΠΟΡΟΟΓ



Patient Informat	tion	
Patient	CORNEAL ANALYZER	G
Patient ID	12345678	E
Date of Birth	01/01/2000	S
	Patient Patient ID	Patient ID 12345678







	Dynamic pupillography				
Dian	Diameter		Pupil Center		
Min	Max	Center Mean	Std Dev.		
3.74 mm	6 26 mm	-0.22 mm	0.02 mm		
3.74 11111	6.36 mm	0.01 mm			
Mesopic pupillography					
Dian	neter	Pupil Center			
Avg. Diam.	Diam. Std Dev.	Center X - Y	Cen. Std Dev.		
4.90 mm	0.33 mm	-0.18 mm	0.33 mm		
4.90 11111	0.55 11111	0.02 mm			
N1 /					
Notes					
I					
CA800 V.1.0.7					

CA800 V.1.0.7		

## **Reports Pupillometry**

# Reports Zernike analysis





CA800 V.1.0.7

# Reports Contact lens fitting

	TOPCON			
Gender Exam Date Surgeon	M 17/02/2015	5 13:3	7:29	
LENSES				
		1		
	Refra	action		
Sphere 125	Cylinder	action Axis 110		VD 12.00
Sphere 1.25	Cylinder -0.50	Axis 110		VD 12.00
1.25	Cylinder -0.50 Lens	Axis 110		12.00
	Cylinder -0.50 Lens Mo	Axis 110		12.00 ase Curve 7.6
1.25 Brand	Cylinder -0.50 Lens Mo Tric	Axis 110 Data odel		12.00 ase Curve



### **Specifications**

Analysed pointsOver 100.000Measured points0 Voer 6.200Corneal coverageUp to 9.8mm on a sphere of radius 8.00mm (42.2 diopters with N=1.3375)Dopter power rangeFrom Ib to 120DResolution4.4'-0.0ID,1 micronAccuracy / Precision axial radius4.4'-0.0IS maltimetric data +/- 2µm at 4mmCapture systemAuto-focus with Auto-captureOutput portsUSB, LANMoltorICD 101 Inch capacitive touch screenPubliometryInternalPubliometryDynamic, Photopic, Mesopic, ScotopicResolutionOro-400C, Relative hundity 30-75% (no dewing), Atmospheric pressure 700-1060/PBPubliometryOro-400C, Relative hundity 30-75% (no dewing), Atmospheric pressure 700-1060/PBPower sourceAC 100-2004 Ar 63 HzPower sourceSolom (M) x 490mm (L), 15 KgConcetionsUSB printer, Network printer, PDF on network shared folder, PDF on USBPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBPrinting optionsSolom (M) x 490mm (L), 15 KgConcetionsUSB printer, Network printer, PDF on network shared folder, PDF on USBPrinting optionsSolom (M) x 490mm (L), 15 KgConcetionsUSB printer, Network printer, PDF on network shared folder, PDF on USBPrinting optionsSolom (M) x 490mm (L), 25 KgConcetionsSolom (M) x 400mm (L), 25 KgPrinting optionsSolom (M) x 400mm (L), 25 KgPower sourceSolom (M) x 400mm (L), 25 KgPower sourceSolom (M) x 400mm (L), 25 KgPower sour	Keratoscope cone	24 rings equally distributed on a 43D sphere
Corneal coverageUp to 9.8mm on a sphere of radius 8.00mm (42.2 diopters with N=1.3375)Diopter power rangeFrom ID to 120DResolution+/- 0.01D, 1 micronAccuracy / Precision axial radius+/- 0.03mm altimetric data +/- 2µm at 4mmCapture systemAuto-focus with Auto-captureOutput portsUSB, LANMonitorLCD 10.1 inch capacitive touch screenDatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFloresceinInage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, ScreenshotPower sourceAC 100-240V 47-63 HzPower consumption320mm (W) x 490mm (H) x 470mm (L), 15 KgPometionsUSB printer, Network printer, PDF on network shared folder, PDF on USBPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBRMZob	Analysed points	Over 100.000
Diopter power rangeFrom Db to 120DResolution.4/- 0.01D, 1 micronAccuracy / Precision axial radius.4/- 0.03mm altimetric data +/- 2µm at 4mmCapture systemAuto-focus with Auto-captureOutput portsUSB, LANMonitorLCD 10.1 inch capacitive touch screenDatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zemike analysis, pupillometry, Toric IOL, ScreenshotVorking environment.100-'40'C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceAC 100-240V 47-63 HzPower consumption.200m (W) x 490mm (H) x 470mm (L), 15 KgConnectionsWi-Fi Optional, LAN integratedPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRM2Gb	Measured points	Over 6.200
Resolution.+/- 0.01D, 1 micronAccuracy / Precision axial radius.+/- 0.03mm altimetric data +/- 2µm at 4mmCapture systemAuto-focus with Auto-captureOutput portsUSB, LANMonitorLCD 10.1 inch capacitive touch screenDatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoKorking environmentOto-Au <sup>o</sup> C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceAC 100-240V 47-63 HzPower consumption<100VAOmnetions320mm (W) x 490mm (H) x 470mm (L), 15 KgConnectionsUSB printer, Network printer, PDF on network shared folder, PDF on USBPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBRMZob	Corneal coverage	Up to 9.8mm on a sphere of radius 8.00mm (42.2 diopters with N=1.3375)
Accuracy / Precision axial radius.+/- 0.03mm altimetric data +/- 2µm at 4mmCapture systemAuto-focus with Auto-captureOutput portsUSB, LANMonitorLCD 10.1 inch capacitive touch screenDatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, ScreenshotVorking environment10°-40°C, Relative humidity 30°75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceAC 100-240V 47-63 HzPomer source320mm (W) x 490mm (H) x 470mm (L), 15 KgConnectionsWi-Fi Optional, LAN integratedPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRM26b	Diopter power range	From 1D to 120D
Capture systemAuto-focus with Auto-captureOutput portsUSB, LANMonitorLCD 10.1 inch capacitive touch screenDatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, ScreenshotPower sourceAC 100-240V 47-63 HzPower consumptionS20mm (W) x 490mm (H) x 470mm (L), 15 KgPinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBPrinting optionsS20bQuestionsWindows 7 EmbeddedRMS20b	Resolution	,+/- 0.01D, 1 micron
Output portsUSB, LANMonitorLCD 101 inch capacitive touch screenDatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, ScreenshotWorking environment10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceAC 100-240V 47-63 HzPower consumption100VADimensions320mm (W) x 490mm (H) x 470mm (L), 15 KgPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRM26b	Accuracy / Precision axial radius	,+/- 0.03mm altimetric data +/- 2µm at 4mm
MonitorLCD 10.1 inch capacitive touch screenDatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, ScreenshotVorking environment10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceAC 100-240V 47-63 HzPower consumption10°UVADimensions320mm (W) x 490mm (H) x 470mm (L), 15 KgPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRAM26b	Capture system	Auto-focus with Auto-capture
DatabaseInternalPupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, ScreenshotWorking environment10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceAC 100-240V 47-63 HzPower consumption<100VADimensions320mm (W) x 490mm (H) x 470mm (L), 15 KgPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRAM2Gb	Output ports	USB, LAN
PupillometryDynamic, Photopic, Mesopic, ScotopicFluoresceinImage, VideoReportCorneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, ScreenshotWorking environment10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceA C 100-240V 47-63 HzPower consumption<100VADimensions320mm (W) x 490mm (L), 15 KgPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRAM26b	Monitor	LCD 10.1 inch capacitive touch screen
Fluorescein       Image, Video         Report       Corneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, Screenshot         Working environment       10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPa         Power source       AC 100-240V 47-63 Hz         Power consumption       <100VA         Dimensions       320mm (W) x 490mm (H) x 470mm (L), 15 Kg         Connections       Wi-Fi Optional, LAN integrated         Printing options       USB printer, Network printer, PDF on network shared folder, PDF on USB         Operating System       Windows 7 Embedded         RAM       2Gb	Database	Internal
Report       Corneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL, Screenshot         Working environment       10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPa         Power source       AC 100-240V 47-63 Hz         Power consumption       <100VA         Dimensions       320mm (W) x 490mm (H) x 470mm (L), 15 Kg         Connections       Wi-Fi Optional, LAN integrated         Power system       USB printer, Network printer, PDF on network shared folder, PDF on USB         Operating System       Windows 7 Embedded         RAM       2Gb	Pupillometry	Dynamic, Photopic, Mesopic, Scotopic
Working environment10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPaPower sourceAC 100-240V 47-63 HzPower consumption<100VADimensions320mm (W) x 490mm (H) x 470mm (L), 15 KgConnectionsWi-Fi Optional, LAN integratedPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRAM2Gb	Fluorescein	Image, Video
Power sourceAC 100-240V 47-63 HzPower consumption<100VA	Report	Corneal map, Comparison map, Contact lens, Height map, Zernike analysis, pupillometry, Toric IOL , Screenshot
Power consumption<100VADimensions320mm (W) x 490mm (L), 15 KgConnectionsWi-Fi Optional, LAN integratedPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRAM2Gb	Working environment	10°-40°C, Relative humidity 30-75% (no dewing), Atmospheric pressure 700-1060hPa
Dimensions     320mm (W) x 490mm (H) x 470mm (L), 15 Kg       Connections     Wi-Fi Optional, LAN integrated       Printing options     USB printer, Network printer, PDF on network shared folder, PDF on USB       Operating System     Windows 7 Embedded       RAM     2Gb	Power source	AC 100-240V 47-63 Hz
ConnectionsWi-Fi Optional, LAN integratedPrinting optionsUSB printer, Network printer, PDF on network shared folder, PDF on USBOperating SystemWindows 7 EmbeddedRAM2Gb	Power consumption	<100VA
Printing options     USB printer, Network printer, PDF on network shared folder, PDF on USB       Operating System     Windows 7 Embedded       RAM     2Gb	Dimensions	320mm (W) x 490mm (H) x 470mm (L), 15 Kg
Operating System     Windows 7 Embedded       RAM     2Gb	Connections	Wi-Fi Optional, LAN integrated
RAM 2Gb	Printing options	USB printer, Network printer, PDF on network shared folder, PDF on USB
	Operating System	Windows 7 Embedded
Hard Disk 250 Gb	RAM	2Gb
	Hard Disk	250 Gb



#### IMPORTANT

Subject to change in design and/or specifications without advanced notice.

#### Topcon Europe Medical B.V.

Essebaan 11; 2908 LJ Capelle a/d IJssel; P.O. Box 145; 2900 AC Capelle a/d IJssel; The Netherlands Phone: +31-(0)10-4585077; Fax: +31-(0)10-4585045 E-mail: medical@topcon.eu; www.topcon-medical.eu

Topcon Danmark Praestemarksvej 25: 4000 Roskilde, Danmark Phone: +45-46-327500; Fax: +45-46-327555 E-mail: info@topcon.dk www.topcon.dk

 Topcon Scandinavia A.B.

 Neongatan 2; P.O. Box 25; 43151 Mölndal, Sweden

 Phone: +46-(0)31-7109200; Fax: +46-(0)31-7109249

 E-mail: medical@topcon.se; www.topcon.se

In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

#### **Topcon España S.A.** HEAD OFFICE; Frederic Mompou, 4;

NEAD OFFICE; Frederic Mompou, 4; 08960 Sant Just Desvern; Barcelona, Spain Phone: +34-93-4734057; Fax: +34-93-4733932 E-mail: medica@topcon.es; www.topcon.es

Topcon Italy Viale dell' Industria 60; 20037 Paderno Dugnano, (MI) Italy Phone: +39-02-9186671; Fax: +39-02-91081091 E-mail: info@topcon.it; www.topcon.it

Topcon France BAT A1; 3 route de la révolte, 93206 Saint Denis Cedex Phone: +33-(0)1-49212323; Fax: +33-(0)1-49212324 E-mail: topcon@topcon.fr; www.topcon-medical.fr Topcon Deutschland GmbH Hanns-Martin-Schleyer Strasse 41; D-47877 Willich, Germany Phone: (+49) 2154-885-0; Fax: (+49) 2154-885-177 E-mail: info@topcon-medical.de; www.topcon-medical.de

**Topcon Polska Sp. z o.o.** ul. Warszawska 23; 42-470 Siewierz; Poland Phone: +48-(0)32-670-50-45; Fax: +48-(0)32-671-34-05 www.topcon-polska.pl

#### Topcon (Great Britain) Ltd.

Topcon House; Kennet Side; Bone Lane; Newbury Berkshire RG14 5PX; United Kingdom Phone: +44-(0)1635-551120; Fax: +44-(0)1635-551170 E-mail: medical@topcon.co.uk, www.topcon.co.uk Topcon Ireland Unit 276, Blanchardstown; Corporate Park 2 Ballycoolin; Dublin 5, Ireland Phone: +353-18975900; Fax: +353-18293915 E-mail: medical@topcon.ie; www.topcon.ie

#### TOPCON TOPCON EUROPE MEDICAL BV 75-1 Hasupuma-cho, Itabashi ku, Tokyo 174-8580, Japa

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, Japan. Phone: 3-3558-2523/2522, Fax: 3-3960-4214, www.topcon.co.jp

con@topcon.fr; www.topcon-medical.fr Ph E-