



EU-1C

Wireless Ophthalmic B scan



BioteK


TECNOLOGIE AVANZATE PER LA DIAGNOSI E LA CHIRURGIA



● Features



Wireless Connect



Easy to be sterilized



Freeze on Probe



Wireless Charging

● SPECIFICATIONS

- Display: Android phone/Tablet.
- Frequency: 8/11 MHz
- Scanning mode: Small Convex
- Scanning depth: 30/40/60/80mm
- Measure mode: Manual
- Angle: 80°@25Hz/112°@12Hz

- Lateral resolution <1mm
- Vertical resolution <0.5mm
- Connection: wifi/USB
- Charging: Wifi/USB
- Battery last: 1.5hrs/Working with charging
- Weight: <150g