



# EXPORT PRICE LIST 2015

● ● ● Medical AG  
Luchten 1262  
CH-9427 Wolfhalden

Tel. +41 71 727 10 50  
Fax +41 71 727 10 55  
info@medical.com  
www.medical.com

**medical**  
SWISS TECHNOLOGY FOR SURGERY

# CONTENT

<b>LENS INJECTION SYSTEM</b>			<b>PAGE</b>
ACCUJECT™ Injector Sets	sterile	single-use	3
ACCUJECT™ SCREW Injector Sets	sterile	single-use	4
NAVIJECT™ Injector Sets	sterile	single-use	5
NAVIJECT™ Cartridge Sets	sterile	single-use	6
NAVIJECT™ Injectors	non-sterile	reusable	6
VISCOJECT™- <b>BIO</b> Injector Sets	sterile	single-use	7
VISCOJECT™- <b>BIO</b> Cartridge Sets	sterile	single-use	8
VISCOJECT™ Injectors	non-sterile	reusable	8
VISCOJECT™ <b>eco</b>	sterile	single-use	9
MULTIJECT™ Cartridges	sterile	single-use	10
MULTIJECT™ Injectors	non-sterile	reusable	10
<b>I/A INSTRUMENTS</b>			
Swiss Mono-manual Systems	non-sterile	reusable	11
Swiss Bi-manual Systems	non-sterile	reusable	11
Swiss Bi-manual Systems	sterile	single-use	12
Swiss Bi-manual Systems	sterile	multiple-use (10 times)	13
Swiss Mono-manual Systems	sterile	single-use	14
Swiss Mono-manual Systems	sterile	multiple-use (10 times)	15
<b>PHACO ACCESSORIES</b>			
Phaco Accessories Sets	sterile	single-use	16
Phaco Tips incl. Sleeves	sterile	multiple-use (10 times)	17
Phaco Tips	non-sterile	reusable	18
COOLFLOW Infusion Sleeves	non-sterile	reusable	19
COOLFLOW Infusion Sleeves	sterile	single-use	20
COMPATIBLES(Cassettes)	sterile	single-use	21
<b>VITREO-RETINAL</b>			
SWISS TROCAR SYSTEMS	sterile	single-use	22
<b>DMEK</b>			
ENDOJECT™	sterile	single-use	23

# ACCUJECT™

## SINGLE-USE INJECTOR SETS, STERILE

consisting of:  
ACCUJECT™ Injection system  
with integrated cartridge\*



for 1-piece IOLs



for 3-piece IOLs

### TOPLOADED - FOR HYDROPHOBIC AND HYDROPHILIC LENSES

BOX OF 10

LP604570	ACCUJECT™ 3.0-1P for 1-piece lenses, inc. size sub 3.2 mm	CHF 120.00
LP604500	ACCUJECT™ 2.6-1P for 1-piece lenses, inc. size sub 2.8 mm	CHF 130.00
LP604530	ACCUJECT™ 2.2-1P for 1-piece lenses, inc. size sub 2.5 mm	CHF 130.00
LP604590	ACCUJECT™ 2.1-1P for 1-piece lenses, inc. size sub 2.4 mm (oval cartridge tip design)	CHF 130.00
LP604510	ACCUJECT™ 2.0-1P for 1-piece lenses, inc. size sub 2.4 mm (round cartridge tip design)	CHF 130.00
LP604520	ACCUJECT™ 2.5-3P for 3-piece lenses, inc. size sub 2.8 mm	CHF 140.00
<b>Designed for MICRO-INCISION lenses only:</b>		
LP604540	ACCUJECT™ 1.8-1P for 1-piece lenses, inc. size sub 2.2 mm	CHF 140.00

### BACKLOADED - FOR HYDROPHOBIC LENSES

BOX OF 10

LP604575	ACCUJECT™ 3.0-BL for 1-piece lenses, inc. size sub 3.2 mm	CHF 120.00
LP604505	ACCUJECT™ 2.6-BL for 1-piece lenses, inc. size sub 2.8 mm	CHF 130.00
LP604535	ACCUJECT™ 2.2-BL for 1-piece lenses, inc. size sub 2.6 mm	CHF 130.00
LP604595	ACCUJECT™ 2.1-BL for 1-piece lenses, inc. size sub 2.4 mm (oval cartridge tip design)	CHF 130.00
LP604515	ACCUJECT™ 2.0-BL for 1-piece lenses, inc. size sub 2.4 mm (round cartridge tip design)	CHF 130.00

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections.  
Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.

#### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	10 %
100	-	199 boxes	15 %
200	-	499 boxes	20 %
500 boxes		and more	25 %

# ACCUJECT™ SCREW

## SINGLE-USE INJECTOR SETS, STERILE

consisting of:  
ACCUJECT™ SCREW Injection system  
with integrated cartridge\*



### BACKLOADED - FOR HYDROPHOBIC LENSES

BOX OF 10

ASB3000	ACCUJECT™ SCREW 3.0-BL for 1-piece lenses, inc. size sub 3.2 mm	CHF 120.00
ASB2600	ACCUJECT™ SCREW 2.6-BL for 1-piece lenses, inc. size sub 2.8 mm	CHF 130.00
ASB2200	ACCUJECT™ SCREW 2.2-BL for 1-piece lenses, inc. size sub 2.5 mm	CHF 130.00
ASB2100	ACCUJECT™ SCREW 2.1-BL for 1-piece lenses, inc. size sub 2.4 mm (oval cartridge tip design)	CHF 130.00
ASB2000	ACCUJECT™ SCREW 2.0-BL for 1-piece lenses, inc. size sub 2.4 mm (round cartridge tip design)	CHF 130.00
<b>Designed for MICRO-INCISION lenses only:</b>		
ASB1800	ACCUJECT™ SCREW 1.8-BL for 1-piece lenses, inc. size sub 2.2 mm	CHF 140.00

### TOPLOADED - FOR HYDROPHOBIC AND HYDROPHILIC LENSES

BOX OF 10

AST2100	ACCUJECT™ SCREW 2.1-TL for 1-piece lenses, inc. size sub 2.4 mm (oval cartridge tip design)	CHF 130.00
AST2000	ACCUJECT™ SCREW 2.0-TL for 1-piece lenses, inc. size sub 2.4 mm (round cartridge tip design)	CHF 130.00
<b>Designed for MICRO-INCISION lenses only:</b>		
AST1800	ACCUJECT™ SCREW 1.8-TL for 1-piece lenses, inc. size sub 2.2 mm	CHF 140.00

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections.  
Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.

#### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	10 %
100	-	199 boxes	15 %
200	-	499 boxes	20 %
500 boxes		and more	25 %

# NAVIJECT™

## SINGLE-USE INJECTOR SETS, STERILE

consisting of: NAVIJECT™ Injector with silicone cushion  
NAVIJECT™ Cartridge\*



<b>WITHOUT LOADING UNIT - FOR HYDROPHOBIC AND HYDROPHILIC LENSES</b>		<b>BOX OF 10</b>
LP604435W	NAVIJECT™ 2.9-1P for 1-piece lenses, inc. size sub 3.2 mm	CHF 105.00
LP604410W	NAVIJECT™ 2.8-1P for 1-piece lenses, inc. size sub 3.0 mm	CHF 105.00
LP604425W	NAVIJECT™ 2.2-1P for 1-piece lenses, inc. size sub 2.5 mm (large loading chamber)	CHF 115.00
LP604426W	NAVIJECT™ 2.2-1P/N for 1-piece lenses, inc. size sub 2.5 mm (small loading chamber)	CHF 115.00

consisting of: NAVIJECT™ Injector with silicone cushion  
NAVIJECT™ Cartridge\*  
NAVIJECT™ Loading unit

<b>WITH LOADING UNIT - FOR HYDROPHOBIC AND HYDROPHILIC LENSES</b>		<b>BOX OF 10</b>
LP604405	NAVIJECT™ 2.9-3P for 3-piece lenses, inc. size sub 3.2 mm with haptic retraction hook	CHF 130.00
LP604435	NAVIJECT™ 2.9-1P for 1-piece lenses, inc. size sub 3.2 mm	CHF 120.00
LP604410	NAVIJECT™ 2.8-1P for 1-piece lenses, inc. size sub 3.0 mm	CHF 120.00
LP604415	NAVIJECT™ 2.7-1P for 1-piece lenses, inc. size sub 2.8 mm	CHF 130.00
LP604425	NAVIJECT™ 2.2-1P for 1-piece lenses, inc. size sub 2.5 mm (large loading chamber)	CHF 130.00
LP604426	NAVIJECT™ 2.2-1P/N for 1-piece lenses, inc. size sub 2.5 mm (small loading chamber)	CHF 130.00

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections.  
Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.

### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	10 %
100	-	199 boxes	15 %
200	-	499 boxes	20 %
500 boxes		and more	25 %

# NAVIJECT™

## SINGLE-USE CARTRIDGE SETS, STERILE

consisting of: NAVIJECT™ Cartridge\*  
NAVIJECT™ Loading unit



### FOR HYDROPHOBIC AND HYDROPHILIC LENSES

BOX OF 10

LC604310	NAVIJECT™ 2.8-1P for 1-piece lenses, inc. size sub 3.0 mm	CHF 85.00
LC604325	NAVIJECT™ 2.2-1P for 1-piece lenses, inc. size sub 2.5 mm	CHF 90.00

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections. Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.

#### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	10 %
100	-	199 boxes	15 %
200	-	499 boxes	20 %
500 boxes		and more	25 %

## REUSABLE INJECTORS, NON-STERILE

### FOR USE WITH NAVIJECT™ CARTRIDGE SETS

1 UNIT

LI604305	NAVIJECT™ Injector, Syringe type „ATRAUMATIC“	CHF 350.00
LI604315	NAVIJECT™ Injector, Screw type „ATRAUMATIC“	CHF 495.00

Syringe type



Screw type



#### DISCOUNT per product line

1	-	9 units	list price
10	-	19 units	10 %
20	-	49 units	20 %
50 units		and more	25 %

# VISCOJECT™-BIO

## SINGLE-USE INJECTOR SETS, STERILE

consisting of: VISCOJECT™ Injector with silicone cushion  
VISCOJECT™-BIO Cartridge\*



### FOR HYDROPHOBIC AND HYDROPHILIC LENSES

BOX OF 10

LP604335C	VISCOJECT™-BIO 2.7 for 1-piece lenses, inc. size sub 2.8 mm	CHF 112.00
LP604340C	VISCOJECT™-BIO 2.2 for 1-piece lenses, inc. size sub 2.5 mm	CHF 112.00
	<b>Designed for MICRO-INCISION lenses only:</b>	
LP604350C	VISCOJECT™-BIO 1.8 for 1-piece lenses, inc. size sub 2.2 mm	CHF 122.00
	<b>Designed for SUB-MICRO-INCISION lenses only:</b>	
LP604360C	VISCOJECT™-BIO 1.5 for 1-piece lenses, inc. size sub 1.8 mm	CHF 130.00

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections.  
Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.

#### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	10 %
100	-	199 boxes	15 %
200	-	499 boxes	20 %
500 boxes		and more	25 %

# VISCOJECT™-BIO

## SINGLE-USE CARTRIDGE SETS, STERILE

consisting of: VISCOJECT™-BIO Cartridge\*  
SILICONE CUSHION



### FOR HYDROPHOBIC AND HYDROPHILIC LENSES

BOX OF 10

LP604235C	VISCOJECT™-BIO 2.7 for 1-piece lenses, inc. size sub 2.8 mm	CHF	82.00
LP604240C	VISCOJECT™-BIO 2.2 for 1-piece lenses, inc. size sub 2.5 mm	CHF	82.00
	<b>Designed for MICRO-INCISION lenses only:</b>		
LP604250C	VISCOJECT™-BIO 1.8 for 1-piece lenses, inc. size sub 2.0 mm	CHF	92.00

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections. Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.

#### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	10 %
100	-	199 boxes	15 %
200	-	499 boxes	20 %
500 boxes		and more	25 %

## REUSABLE INJECTORS, NON-STERILE

### FOR USE WITH VISCOJECT™-BIO CARTRIDGE SETS

1 UNIT

LI604205	VISCOJECT™ Injector, Syringe type „ATRAUMATIC“	CHF	425.00
LI604215	VISCOJECT™ Injector, Screw type „ATRAUMATIC“	CHF	495.00

Syringe type



Screw type



#### DISCOUNT per product line

1	-	9 units	list price
10	-	19 units	10 %
20	-	49 units	20 %
50 units		and more	25 %



# VISCOJECT™ **eco**



SINGLE-USE INJECTOR-SETS, STERILE

consisting of: VISCOJECT™ **eco** Injector with silicone cushion  
VISCOJECT™ Cartridge\*

**FOR HYDROPHOBIC AND HYDROPHILIC LENSES**

BOX OF 20

VE2700	VISCOJECT™ <b>eco</b> 2.7 for 1-piece lenses, inc. size sub 2.8 mm	CHF 120.00
VE2200	VISCOJECT™ <b>eco</b> 2.2 for 1-piece lenses, inc. size sub 2.5 mm	CHF 120.00

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections.  
Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.



**DISCOUNT per product line**

1	-	49 boxes	list price
50	-	99 boxes	10 %
100	-	199 boxes	15 %
200 boxes		and more	20 %

# MULTIJECT™

## SINGLE-USE CARTRIDGES, STERILE

consisting of: MULTIGLIDE™ Cartridge\*



### FOR HYDROPHOBIC AND HYDROPHILIC LENSES

1 UNIT

LC604210	OVALGLIDE™ Cartridge for 1- and 3-piece lenses, inc. size sub 3.2 mm	CHF	4.90
LC604215	MACROGLIDE™ Cartridge for 1-piece lenses, inc. size sub 3.5 mm	CHF	4.90
LC604220M	MICROGLIDE™ Cartridge for 1- and 3-piece lenses, inc. size sub 3.0 mm	CHF	5.60
LS604500	SILICONE CUSHION for multiple-use injectors LI604205 and LI604215	CHF	3.50

\* Incision sizes depend on surgical techniques. Listed incision sizes are for „into-the bag“-injections. Injection techniques „into-the-wound“ or „wound-assisted“ result in smaller incisions.

#### DISCOUNT per product line

1	-	199 units	list price
200	-	499 units	10 %
500	-	1'999 units	20 %
2000		and more	25 %

## REUSABLE INJECTORS, NON-STERILE

### FOR USE WITH MULTIJECT™ CARTRIDGES

1 UNIT

LI604200	MULTIJECT™ Injector, Syringe-Type „GROOVE“ for 1- and 3-piece lenses	CHF	425.00
LI604210	MULTIJECT™ Injector, Screw-Type „GROOVE“ for 1- and 3-piece lenses	CHF	495.00
LI604205	MULTIJECT™ Injector, Syringe-Type „ATRAUMATIC“ for 1-piece lenses for use with LS604500 silicone cushion	CHF	425.00
LI604215	MULTIJECT™ Injector, Screw-Type „ATRAUMATIC“ for 1-piece lenses for use with LS604500 silicone cushion	CHF	495.00



#### DISCOUNT per product line










1	-	9 units	list price
10	-	19 units	10 %
20	-	49 units	20 %
50 units		and more	25 %

# I/A INSTRUMENTS

## REUSABLE SWISS MONO- AND BI-MANUAL I/A SYSTEMS


### SWISS QUICK SWITCH MONO-MANUAL I/A INSTRUMENTS

1 UNIT

SQS100	Lightweight I/A handpiece		CHF 235.00
SQS110	I/A tip, straight Ø 1.80 mm, aspiration port 0.35 mm		CHF 180.00
SQS115	I/A polished tip, straight, Ø 1.40 mm, asp. port 0.35 mm		CHF 190.00
SQS120	I/A tip, 45° angled, Ø 1.80 mm, asp. port 0.35 mm		CHF 198.00
SQS130	I/A tip, 90° angled, Ø 1.80 mm, asp. port 0.35 mm		CHF 198.00
SQS140	I/A tip, curved, Ø 1.80 mm, aspiration port 0.35 mm		CHF 198.00
SQS150	I/A tip, Simcoe Style, gold-plated, asp. port 0.35 mm		CHF 250.00
SQS116	I/A tip, straight, 19G, asp. port 0.35 mm with silicone sleeve		CHF 180.00
SQS119	I/A tip, 45° angled, 19G, asp. port 0.35 mm with silicone sleeve		CHF 190.00


### SWISS BI-MANUAL I/A ASPIRATION INSTRUMENTS

1 UNIT

SBA100A	Asp. handpiece <b>21G</b> , aeroplane shape, with polished tip, asp. port 0.35 mm		CHF 140.00
SBA103A	Asp. handpiece <b>21G</b> , aeroplane shape, with polished tip, asp. port 0.40 mm		CHF 140.00
SBA104A	Asp. handpiece <b>21G</b> , aeroplane shape, with textured tip, asp. port 0.40 mm		CHF 160.00
SBA105A	Asp. handpiece <b>21G</b> , aeroplane shape, with textured tip, asp. port 0.35 mm		CHF 160.00
SBA101B	Asp. handpiece <b>23G</b> , bullet shape, with polished tip, asp. port 0.30 mm		CHF 190.00
SBA106B	Asp. handpiece <b>23G</b> , bullet shape, with textured tip, asp. port 0.30 mm	CHF 210.00	
SBA110	Asp. handpiece, with luer adapter for aspiration cannulas	CHF 110.00	

### SWISS BI-MANUAL I/A IRRIGATION INSTRUMENTS

1 UNIT

SBI120A	Irr. handpiece <b>21G</b> , aeroplane shape, with polished tip, irr. port 0.50 mm		CHF 140.00
SBI125A	Irr. handpiece <b>21G</b> , aeroplane shape, with textured tip, irr. port 0.50 mm		CHF 160.00
SBI121B	Irr. handpiece <b>23G</b> , bullet shape, with polished tip, irr. port 0.30 mm	CHF 190.00	
SBI126B	Irr. handpiece <b>23G</b> , bullet shape, with textured tip, irr. port 0.30 mm	CHF 210.00	
SBI130	Irr. handpiece, with luer adapter for irrigation cannulas	CHF 110.00	

# I/A INSTRUMENTS


SINGLE-USE SWISS BI-MANUAL I/A SYSTEMS, STERILE




## SWISS BI-MANUAL I/A SYSTEMS FOR THE BI-MANUAL TECHNIQUE

BOX OF 10

SBS105	<p>Aspiration-handpiece 0.8 mm (21G) with curved cannula, aspiration port 0.35 mm, textured tip for capsule polishing AEROPLANE nose shape</p> <p>Irrigation-handpiece 0.8 mm (21G) with curved cannula, irrigation port 0.50 mm for active rinsing</p> <p>For incisions of 1.1 mm</p>	CHF 110.00
SBS110A	<p>Aspiration-handpiece 0.8 mm (21G) with curved cannula, aspiration port 0.35 mm, textured tip for capsule polishing, AEROPLANE nose shape</p> <p>Irrigation-handpiece 0.8 mm (21G) with curved cannula, two irrigation ports 0.45 mm, AEROPLANE nose shape</p> <p>For incisions of 1.1 mm</p>	CHF 140.00



SBS105



SBS110A

### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	5 %
100	-	199 boxes	10 %
200 boxes		and more	15 %

# I/A INSTRUMENTS

## MULTIPLE-USE SWISS BI-MANUAL I/A SYSTEMS, STERILE


Sets delivered sterile, components can be autoclaved at 121°C/20 min or 134°C/5 min up to 10 times.




### SWISS BI-MANUAL I/A SYSTEMS

BOX OF 10

<p><b>SBS105RU</b></p> <p>Aspiration-handpiece 0.8 mm (21G) with curved cannula, aspiration port 0.35 mm, with textured aspiration head for capsule polishing AEROPLANE nose shape</p> <p>Irrigation-handpiece 0.8 mm (21G) with curved cannula, two irrigation ports 0.5 mm for active rinsing</p> <p>For incisions of 1.1 mm</p>	<p>CHF 160.00</p>
<p><b>SBS110ARU</b></p> <p>Aspiration-handpiece 0.8 mm (21G) with curved cannula, aspiration port 0.35 mm, with textured aspiration head for capsule polishing AEROPLANE nose shape</p> <p>Irrigation-handpiece 0.8 mm (21G) with curved cannula, two irrigation ports 0.45 mm AEROPLANE nose shape</p> <p>For incisions of 1.1 mm</p>	<p>CHF 190.00</p>



SBS105RU



SBS110ARU

#### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	5 %
100	-	199 boxes	10 %
200 boxes		and more	15 %

# I/A INSTRUMENTS


SINGLE-USE SWISS MONO-MANUAL I/A SYSTEMS, STERILE



## SWISS MONO-MANUAL I/A SYSTEMS FOR THE CO-AXIAL TECHNIQUE

BOX OF 10

SMS120	I/A handpiece with aspiration cannula <b>angled 45° / 19G</b> , for incisions of 2.8 mm, aspiration port 0.35 mm, with textured aspiration head for capsule polishing	CHF 120.00
SMS130	I/A handpiece with aspiration cannula <b>angled 65° / 19G</b> , for incisions of 2.8 mm, aspiration port 0.35mm, with textured aspiration head for capsule polishing	CHF 120.00
SMS115	I/A handpiece with aspiration cannula <b>straight / 20G</b> , for incisions of 2.5 mm, aspiration port 0.35 mm, with textured aspiration head for capsule polishing	CHF 100.00
SMS125P	I/A handpiece with aspiration cannula <b>angled 45° / 20G</b> , for incisions of 2.5 mm, aspiration port 0.35 mm, with polished aspiration head	CHF 120.00
SMS135	I/A handpiece with aspiration cannula <b>angled 65° / 20G</b> , for incisions of 2.5 mm, aspiration port 0.35 mm, with textured aspiration head for capsule polishing	CHF 120.00
SMS155P	I/A handpiece with aspiration cannula <b>angled 45° / 21G</b> , for incisions of 2.2 mm, aspiration port 0.35 mm, with polished aspiration head	CHF 130.00
SMS160P	I/A handpiece with aspiration cannula <b>curved / 21G</b> , for incisions of 2.2 mm, aspiration port 0.35 mm, with polished aspiration head	CHF 130.00
SMS170P	I/A handpiece with aspiration cannula <b>angled 45° / 23G</b> , for incisions of 2.0 mm, aspiration port 0.30 mm, with polished aspiration head	CHF 140.00
SMS175P	I/A handpiece with aspiration cannula <b>curved / 23G</b> , for incisions of 2.0 mm, aspiration port 0.30 mm, with polished aspiration head	CHF 140.00



### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	5 %
100	-	199 boxes	10 %
200 boxes		and more	15 %

# I/A INSTRUMENTS

## MULTIPLE-USE SWISS MONO-MANUAL I/A SYSTEMS, STERILE

Sets delivered sterile, components can be autoclaved at 121°C/20 min or 134°C/5 min up to 10 times.



SWISS MONO-MANUAL I/A SYSTEMS FOR THE CO-AXIAL TECHNIQUE		BOX OF 10
SMS120RU	I/A handpiece with aspiration cannula <b>angled 45° / 19G</b> asp. port 0.35 mm, with textured aspiration head for capsule polishing for incisions of 2.8 mm	CHF 170.00
SMS125PRU	I/A handpiece with aspiration cannula <b>angled 45° / 20G</b> asp. port 0.35 mm, with polished aspiration head for incisions of 2.5 mm	CHF 170.00
SMS135RU	I/A handpiece with aspiration cannula <b>angled 65° / 20G</b> asp. port 0.35 mm, with textured aspiration head for capsule polishing for incisions of 2.5 mm	CHF 170.00
SMS170PRU	I/A handpiece with aspiration cannula <b>angled 45° / 23G</b> asp. port 0.30 mm, with polished aspiration head for incisions of 2.0 mm	CHF 190.00

### DISCOUNT per product line

1	-	49 boxes	list price
50	-	99 boxes	5 %
100	-	199 boxes	10 %
200 boxes		and more	15 %

# PHACO-ACCESSORIES

## SINGLE-USE PHACO-ACCESSORIES SETS, STERILE











consisting of: 1 Titanium Tip with tip-wrench  
1 Infusion Sleeve  
1 Test Chamber



compatible with

5R + 4K-Hub: B+L/Storz, AMO, Alcon, Nidek, Oertli, Ruck, Geuder, D.O.R.C. and others

XXK-Hub: Alcon Legacy® Series 20000, Accurus®, Infinity®

					
		5R-Hub	4K-Hub	XXK-Hub	BOX OF 10
	CLASSIC 19G, 30° for 2.8 mm incisions	PA603000	PA603500	PA603300	CHF 225.00
	MICROSONIC 20G, 30° for 2.5 mm incisions	PA603100	PA603510	PA603310	CHF 260.00
	HYPERSONIC 19/20G, 30° for 2.8 mm incisions	PA603200	PA603520	PA603350	CHF 260.00
	NANO 20/21G, 30° (flared) for 2.2 mm incisions	PA603190	PA603530	PA603320	CHF 285.00
	SUBSONIC 21G, 30° for 2.2 mm incisions	PA603070	PA603540	PA603390	CHF 285.00
	COMICS-SUB2 23G, 30° for 1.8 mm incisions	PA603020	PA603550		CHF 380.00
					1 UNIT
	TIP-WRENCH (non-sterile)	PK602402	PK602401	PK602401	CHF 65.00

Upon request, phaco tips are also available in 0°, 15°, 45° and/or curved version and can be completed with a bimanual or monomanual I/A system.  
Minimum order quantity for these products is 10 boxes per order.

### DISCOUNT per product line

1	-	9 boxes	list price
10	-	49 boxes	10 %
50	-	and more	20 %



# PHACO-ACCESSORIES

## MULTIPLE-USE PHACO TIPS AND SLEEVES, STERILE

consisting of: 1 Titanium Tip with autoclavable tip-wrench  
1 Infusion Sleeve













compatible with

5R + 4K-Hub: B+L/Storz, AMO, Alcon, Nidek, Oertli, Ruck, Geuder, D.O.R.C. and others

XXK-Hub: Alcon Legacy® Series 20000, Accurus®, Infinity®

Sets delivered sterile, components can be autoclaved at 121°C/20 min up to 10 times.

				
	5R-Hub	4K-Hub	XXK-Hub	BOX OF 10
 CLASSIC 19G, 30° for 2.8 mm incisions	PA603000RU	PA603500RU	PA603300RU	CHF 220.00
 MICROSONIC 20G, 30° for 2.5 mm incisions	PA603100RU	PA603510RU	PA603310RU	CHF 257.00
 HYPERSONIC 19/20G, 30° for 2.8 mm incisions	PA603200RU	PA603520RU	PA603350RU	CHF 257.00
 NANO 20/21G, 30° (flared) for 2.2 mm incisions	PA603190RU	PA603530RU	PA603320RU	CHF 280.00
 SUBSONIC 21G, 30° for 2.2 mm incisions	PA603070RU	PA603540RU	PA603390RU	CHF 280.00
 COMICS-SUB2 23G, 30° for 1.8 mm incisions	PA603020RU	PA603550RU		CHF 360.00
				1 UNIT
 TIP-WRENCH (non-sterile)	PK602402	PK602401	PK602401	CHF 65.00

Upon request, phaco tips are also available in 0°, 15°, 45° and/or curved version.  
Minimum order quantity for these products is 10 boxes per order.

### DISCOUNT per product line

1	-	9 boxes	list price
10	-	49 boxes	10 %
50	-	and more	20 %

# PHACO-ACCESSORIES

## REUSABLE PHACO TIPS, NON-STERILE











consisting of: 6 Titanium Tips

compatible with

5R + 4K-Hub: B+L/Storz, AMO, Alcon, Nidek, Oertli, Ruck, Geuder, D.O.R.C. and others

XXK-Hub: Alcon Legacy® Series 20000, Accurus®, Infinity®



					
		5R-Hub	4K-Hub	XXK-Hub	BOX OF 6
	CLASSIC 19G, 30° for 2.8 mm incisions	PT602200B	PT602500B	PT602300B	CHF 112.00
	MICROSONIC 20G, 30° for 2.5 mm incisions	PT602210B	PT602510B	PT602310B	CHF 144.00
	HYPERSONIC 19/20G, 30° for 2.8 mm incisions	PT602250B	PT602520B	PT602350B	CHF 144.00
	NANO 20/21G, 30° (flared) for 2.2 mm incisions	PT602240B	PT602530B	PT602340B	CHF 200.00
	SUBSONIC 21G, 30° for 2.2 mm incisions	PT602070B	PT602540B	PT602390B	CHF 200.00
	COMICS-SUB2 23G, 30° for 1.8 mm incisions	PT602220B	PT602550B		CHF 280.00
					1 UNIT
	TIP-WRENCH (non-sterile)	PK602402	PK602401	PK602401	CHF 65.00

Upon request, phaco-tips are also available in 0°, 15°, 45° and/or curved version.  
Minimum order quantity for these products is 10 boxes per order.

### DISCOUNT per product line

1	-	9 boxes	list price
10	-	49 boxes	10 %
50	-	and more	20 %


# PHACO-ACCESSORIES

## REUSABLE COOLFLOW INFUSION SLEEVES, NON-STERILE



### TO BE USED WITH MEDICEL PHACO TIPS

BAG OF 6

CF602300	COOLFLOW 19G	yellow, for 2.8 mm incisions	CHF	21.00
CF602301	COOLFLOW 19G	blue, for 2.8 mm incisions	CHF	21.00
CF602310	COOLFLOW 20G	beige, for 2.5 mm incisions	CHF	27.00
CF602340	COOLFLOW 21G	clear, for 2.2 mm incisions	CHF	43.00
CF602350	COOLFLOW 23G	orange, for 1.8 mm incisions	CHF	45.00
TC602450	TEST CHAMBER		CHF	19.00


BOX OF 4

CF602320	COOLFLOW <sup>PLUS</sup> 19G	Rigid no-burn, yellow	CHF	48.00
CF602330	COOLFLOW <sup>PLUS</sup> 20G	Rigid no-burn, beige	CHF	52.00

DISCOUNT per product line  
50 bags/boxes and more 10 %

### TO BE USED WITH MEDICEL PHACO TIPS

BAG OF 500

CF602300B	COOLFLOW 19G	yellow, for 2.8 mm incisions	CHF	1'575.00
CF602301B	COOLFLOW 19G	blue, for 2.8 mm incisions	CHF	1'575.00
CF602310B	COOLFLOW 20G	beige, for 2.5 mm incisions	CHF	2'025.00
CF602340B	COOLFLOW 21G	clear, for 2.2 mm incisions	CHF	3'240.00
TC602450B	TEST CHAMBER		CHF	1'440.00

DISCOUNT per product line  
1 - 3 bags list price  
4 bags and more 20 %

# PHACO-ACCESSORIES

SINGLE-USE PHACO-ACCESSORIES WITHOUT PHACO TIPS, STERILE

## INFUSION SLEEVES AND TEST CHAMBER BOX OF 10

PA603901	Accessories Set, 19G (COOLFLOW sleeve, blue and test chamber)	CHF 60.00
PA603902	Accessories Set, 20G small port (COOLFLOW sleeve, beige and test chamber)	CHF 70.00
PA603903	Accessories Set, 21G small port (COOLFLOW sleeve, clear and test chamber)	CHF 100.00

## INFUSION SLEEVES BOX OF 10

PA603951	COOLFLOW 19G Sleeve, blue, for 19G phaco tips	CHF 31.00
PA603952	COOLFLOW 20G Sleeve, beige, for 20G phaco tips	CHF 40.00
PA603953	COOLFLOW 21G Sleeve, clear, for 21G phaco tips	CHF 65.00

## TEST CHAMBER BOX OF 10

TC602451	TEST CHAMBER	CHF 28.00
----------	--------------	-----------



### DISCOUNT per product line

1	-	49 boxes	list price
50	-	199 boxes	10 %
200	-	499 boxes	15 %
500 boxes		and more	20 %

# COMPATIBLES

## SINGLE-USE ULTRA-HIGHVAC XXK CASSETTE SYSTEMS, STERILE

consisting of: 1 XXK ULTRA-HIGHVAC cassette  
1 I/A tubing set  
1 Administration set (with filter)  
1 Tray support cover

compatible with  
Alcon Legacy® Series 20000



BOX OF 12

XXK100	XXK ULTRA-HIGHVAC SET without sleeves	CHF 405.00
XXK101	XXK ULTRA-HIGHVAC SET for 19G phaco tips with COOLFLOW 19G sleeve, blue and test chamber	CHF 424.00
XXK120	with XXK MICROSONIC 20G, 30° phaco tip and 20G sleeve, beige for 2.5 mm incision	CHF 600.00

## SINGLE-USE CONCEPT CASSETTE SYSTEMS, STERILE

compatible with  
Bausch + Lomb/Storz Systems Millennium™, Premiere®, Protégé®, Daisy®



BOX OF 15

CPT100	1 Single-use CONCEPT cassette 1 Single-use I/A Tubing set	CHF 270.00
BOX OF 10		
CPT105	1 Single-use CONCEPT cassette 1 Single-use I/A Tubing set 1 Single-use filtered administration set	1 Single-use infusion sleeve 19G, blue 1 Single-use test chamber 2 Single-use tip-wrenches (5R, 4K) CHF 250.00
CPT109	1 Single-use CONCEPT cassette 1 Single-use I/A Tubing set 1 Single-use filtered administration set 1 Single-use tray support cover	1 Remote control bag 1 Single-use machine cover 1 Single-use infusion sleeve 19G, blue 1 Single-use test chamber CHF 300.00

# VITREO-RETINAL

## SINGLE-USE SWISS TROCAR SYSTEMS, STERILE



### SWISS TROCAR SYSTEMS

BOX OF 10

STS100	<p>SWISS TROCAR SYSTEM, 23G, ONE-STEP WITH PLUGS including</p> <ul style="list-style-type: none"> <li>- 2 Ultra sharp 23G classic blade handpieces with highly polished trocars for low instrument friction</li> <li>- 1 Incisional template</li> <li>- 1 23G blade handpiece with infusion cannula</li> <li>- 1 23G Infusion terminal (BSS) with luer lock connection</li> <li>- 3 Scleral plugs</li> <li>- 1 Holding forceps adjusted to the shape of the trocar</li> <li>- 1 Speculum</li> </ul>	CHF 990.00
STS101	<p>SWISS TROCAR SYSTEM, 23G, ONE-STEP WITH PLUGS including</p> <ul style="list-style-type: none"> <li>- 3 Ultra sharp 23G classic blade handpieces with highly polished titanium trocars for minimal instrument friction</li> <li>- 1 Incisional template</li> <li>- 1 23G Infusion terminal (BSS) with pre-mounted trocar adapter and luer lock connection</li> <li>- 3 Scleral plugs</li> <li>- 1 Holding forceps adjusted to the shape of the trocar</li> </ul>	CHF1090.00
STS102	<p>SWISS TROCAR SYSTEM, 23G, ONE-STEP WITH VALVE including</p> <ul style="list-style-type: none"> <li>- 3 Ultra sharp 23G MVR-blade handpieces with highly polished titanium trocars for minimal instrument friction and self-sealing valve</li> <li>- 1 Incisional template</li> <li>- 1 23G Infusion terminal (BSS) with pre-mounted trocar adapter and luer lock connection</li> <li>- 1 Holding forceps adjusted to the shape of the trocar</li> </ul>	CHF1090.00
STS250	<p>SWISS TROCAR SYSTEM, 25G, ONE-STEP WITH VALVE including</p> <ul style="list-style-type: none"> <li>- 3 Ultra sharp 25G MVR blade handpieces with highly polished titanium trocars for minimal instrument friction and self-sealing valve</li> <li>- 1 Incisional template</li> <li>- 1 25G infusion terminal (BSS) with pre-mounted trocar adapter and luer lock connection</li> <li>- 1 Holding forceps adjusted to the shape of the trocar</li> </ul>	CHF 1130.00
STS260	<p>SWISS TROCAR SYSTEM, 25G, ONE-STEP WITH VALVE including</p> <ul style="list-style-type: none"> <li>- 3 Ultra sharp 25G MVR-blade handpieces with highly polished trocars for minimal instrument friction and self-sealing valve</li> <li>- 1 Incisional template</li> <li>- 1 25G infusion terminal (oil)</li> <li>- 1 Holding forceps adjusted to the shape of the trocar</li> </ul>	CHF 1130.00

**DISCOUNT per product line**

1	-	5 boxes	list price
6	-	19 boxes	10 %
20 boxes		and more	15 %

# ENDOJECT™

SINGLE-USE DMEK INJECTOR SET, STERILE

consisting of: ENDOJECT™ Injector with silicone cushion  
 ENDOJECT™ Cartridge



BOX OF 1

EJ2200	ENDOJECT™ 2.2 for safe implantation of grafts during Descemet Membrane Endothelial Keratoplasty	CHF 49.00
--------	---	-----------

DISCOUNT per product line			
1	-	99 boxes	list price
100	-	199 boxes	10 %
200	-	399 boxes	15 %
400 boxes		and more	20 %