



**BIOTECH**  
OPHTHALMICS

# CATARACT SURGERY SOLUTIONS

[www.biotechvisioncare.com](http://www.biotechvisioncare.com)

## Bio-Hyalur Range Sodium Hyaluronate Ophthalmic Solution

The best options for the protection & control you need in ophthalmic surgery

- Cataract Extraction
- Intraocular Lens (IOL) implantation
- Corneal transplantation surgery
- Glaucoma filtering surgery

### Advantages:

- Source - Bacterial fermentation
- Ideal for lubricating the cartridges of IOLs for implantation
- Excellent endothelial protection
- Easy to inject and remove from the eye
- High viscosity for further expansion of the anterior chamber & higher maintenance of spaces while doing ophthalmic surgery
- Latex free

## Bio-Hyalur

Sodium Hyaluronate 1% Ophthalmic Solution

### Description:


Bio-Hyalur is a highly purified Sodium Hyaluronate 10 mg/ml to use as ophthalmic visco-elastic solution during Anterior Segment surgeries. It possesses viscosity of 300,000 mPas @ Zero shear rate.

### Indications:

- Intraocular lens implantation & Cataract surgery
- Glaucoma surgery
- Corneal transplantation

### Packaging:

- Bio-Hyalur is supplied in 1.0 ml in sterile PFS along with 27 gauge sterile cannula



## Bio-Hyalur Plus

Sodium Hyaluronate 1.4% Ophthalmic Solution

### Description:

Bio-Hyalur Plus is a highly purified Sodium Hyaluronate 14 mg/ml to use as ophthalmic visco-elastic solution during Anterior Segment surgeries. It possesses viscosity of 400,000 mPas @ Zero shear rate.

### Indications:

- Intraocular lens implantation & Cataract surgery
- Glaucoma surgery
- Corneal transplantation

### Packaging:

- Bio-Hyalur Plus is supplied in 1.0 ml in sterile PFS along with 27 gauge sterile cannula



## Bio-Hyalur HV

Sodium Hyaluronate 2.4% Ophthalmic Solution

### Description:

Bio-Hyalur HV is a highly purified Sodium Hyaluronate 24 mg/ml to use as ophthalmic visco-elastic solution during Anterior Segment surgeries. It possesses viscosity of 1,000,000 mPas @ Zero shear rate.

### Indications:

- Anterior segment surgery
- Corneal transplantation
- Phacoemulsification in complicated cataract surgery
- Paediatric cataract surgery

### Packaging:

- Bio-Hyalur HV is supplied in 1.0 ml in sterile PFS along with 27 gauge sterile cannula



## Bio-Hyalur SV

Sodium Hyaluronate 3.0% Ophthalmic Solution

### Description:

Bio-Hyalur SV is a highly purified Sodium Hyaluronate 30 mg/ml to use as ophthalmic visco-elastic solution during Anterior Segment surgeries. It possesses viscosity of 1,300,000 mPas @ Zero shear rate.

### Indications:

- Anterior segment surgery
- Corneal transplantation
- Phacoemulsification in complicated cataract surgery
- Paediatric cataract surgery

### Packaging:

- Bio-Hyalur SV is supplied in 1.0 ml in sterile PFS along with 27 gauge sterile cannula



## Bio-Hyalur DUO

Combi Pack of Bio-Hyalur 1.0% & Bio-Hyalur HV 2.4%

- Bio-Hyalur DUO offer unique properties that work together and give the ultimate protection, control, clarity, and ease-of-use for every cataract and phacoemulsification procedure.
- The combination allows for ideal viscoelastic utilization during each phase of surgery which the surgeon looks for.

### Packaging:

We offer Bio-Hyalur 1.0% 1.0ml & Bio-Hyalur HV 2.4% 1.0ml in a single pack.



## CATARACT SURGERY SOLUTIONS



### Factory

1st Floor, 555 / 556 / 557 , Nr. Shubham Tex-O-Pack,  
Khatraj, Dist. Gandhinagar, Gujarat, India.

### India Business

210, Sarthik - II, Opp. Rajpath Club,  
Sarkhej-Gandhinagar Highway,  
Ahmedabad-380054, Gujarat, India.  
**Ph. / Fax:** 079-40038374,32420133  
**E-mail:** sales@biotechvisioncare.com

### International Business

401, Sarthik- II, Opp. Rajpath Club,  
Sarkhej-Gandhinagar Highway,  
Ahmedabad-380054, Gujarat, India.  
**Ph.:** +91-79-26870896/97/98 **Fax:** +91-79-40033767  
**E-mail:** intbusiness@biotechvisioncare.com